

Model Number

K9907C

# MEDIUM PRESSURE CALIBRATION SYSTEM

Revision: NR

ECN#:

## GENERAL

Excitation Range	100 to 1000 psi [70 to 700 kPa]
Calibration Technique	Dynamic Step Pulse
Typical Rise Time	20 to 40 microseconds
Measurements Supported	Sensitivity, Amplitude Linearity
Pressure Units Supported	psi, Pa, kPa, Mpa
Sensor Types Supported	Dynamic Pressure
Operation Types Supported	IEPE, Charge, Voltage
Operation Types Supported (Requires Optional Hardware)	Piezoresistive (Bridge), Capacitive, Differential Charge
Included PC	Windows 10
Test Parameters	Managed by system software
Test Pass/Fail Criteria	Managed by system software
Data Acceptance	Manual or automatic by system software
Data Management	Via User Interface or ODBC SQL
Report Generation	Automatic with MS Excel

## DATA ACQUISITION SYSTEM

Resolution	24 Bits
ADC Type	Delta-Sigma

## INCLUDED SUT SIGNAL CONDITIONER – ICP®/CHARGE


Manufacturer	PCB
Model	443B02
Input Type	ICP/Charge
Input Connector	BNC Jack
Output Type	Voltage <sup>[1]</sup>
Output Connector	BNC Jack
Variable Gain Increment	0.1 to 1000
Output Range	± 10 Volts
Control Interface	Front Panel, RS-232

## MISCELLANEOUS

The K9907C integrates software and hardware and is statistically verified before shipment and after arriving at final destination.

## MAJOR COMPONENTS

Aronson Dynamic Pressure Calibrator, Windows PC, ICP/Charge Signal Conditioning, Data Acquisition Hardware, Software, and Verification Sensor



NOTES:  
[1] Output Impedance <1 Ohm

All specifications are at room temperature unless otherwise specified.

ICP® is a registered trademark of PCB Piezotronics, Inc.

In the interest of constant product improvement, specifications may change without notice.

Project Engineer:

ATG

Product Manager:

IND

SAM Team Leader:

CAO

Spec Number:

PS-0162

Date: 11/29/2020

Date: 11/24/2020

Date: 12/1/20

Page 1 of 1

**THE MODAL SHOP**  
AN AMPHENOL COMPANY

10310 Aerohub Boulevard  
Cincinnati, OH 45215, USA

info@modalshop.com  
+1 513.351.9919

+1 800.860.4867  
Fax: +1 513.458.2172

SAM-F020 revC 10/20/20